



PacificDx
Seamless Results.

Novel four color FISH probe for simultaneous detection of RET break-apart and PDGFR α amplification

Aditi Khurana , Suman Verma, Philip D. Cotter, and Mathew W. Moore
PacificDx, 5 Mason, Irvine, CA

INTRODUCTION

A number of critical genomic rearrangements are associated with therapeutics in non small cell lung carcinoma (NSCLC). There is frequently limited FFPE material available and in the interest of utilizing scarce clinical samples, we developed a four-color RET/PDGFR α /KIT/CC4 FISH probe. The probe provides coverage of the PDGFR α gene (4q12), the KIT gene (4q12), a copy control probe located at 4p12 on chromosome 4 and detection of RET gene rearrangements (10q11.2). This novel four color FISH probe has been validated for clinical use on 4 to 5 um thick FFPE sections.

CONCLUSION

An assay validation was performed using 23 negative samples consisting of de identified Human Bio specimens, xenografts and cell lines for genes of interest including 2 positive samples for PDGFR α amplification and 2 for RET rearrangement. The assay demonstrated 100% concordance with known standards and yielded clinically actionable cut-offs that comply with industry standards.

PROBE DESIGN




Figure 1 : Probe design for the RET/PDGFR α FISH probe manufactured by Biocare Medical LLC for RUO purposes only. The RET proto-oncogene encodes a receptor tyrosine kinase, and chromosomal rearrangements that generate a fusion gene resulting in the juxtaposition of the C-terminal region of the RET protein with an N-terminal portion of another protein, can also lead to constitutive activation of the RET kinase. The RET gene is localized to chromosome 10 (10q11.2) and has 21 exons. PDGFR α or the Platelet- derived growth factor receptor mutations are associated with idiopathic hypereosinophilic syndrome, somatic, familial gastrointestinal stromal tumors, and a variety of other cancers. The PDGFR α gene is localized to chromosome 4 (4q12) and has 28 exons. This gene plays a role in organ development, wound healing, and tumor progression.

PROBE LOCALIZATION




Figure 2: Two known normal male metaphase slides were prepared using standard techniques and hybridized using standard methodologies. Above is the fluorescent image seen under a Zeiss microscope compared to an inverted DAPI image that allows identification of chromosomes as a quality control measure. A total of 50 metaphases were analyzed, and five representative images were captured to confirm proper probe localization. The PDGFR α /KIT/CC4/RET novel four color probe set displayed 100% probe sensitivity and 100% probe specificity and there was no evidence of background or cross hybridization detected.

ACCURACY

Sample ID/ Name	Results from Validation	
	Result by FISH	Result by NGS
NCIH 1693	5.62 copies PDGFR α Negative for RET rearrangement	4.6 copies PDGFR α RET mutation Not Detected
LC2Ad	Negative for PDGFR α amplification RET rearrangement >68%	PDGFR α amplification Not Detected CCD6-RET translocation Detected
GM19204 (Hapmap, WT)	Negative for PDGFR α amplification Negative for RET gene rearrangement	PDGFR α amplification Not Detected RET mutation Not Detected

Figure 3: Accuracy of FISH assay: all samples were assayed according to SOP approved by the Medical Director. For this validation, the normal cohort consisted of sectioned FFPE cell-line derived blocks and NSCLC specimens previously screened using an orthogonal technology at PacificDx. The abnormal cohort comprised of known positive NSCLC samples, sectioned FFPE cell-line derived blocks, and custom xenografts . Assay acceptance criteria: >95% concordance in between Negative/ Positive calls between NGS and FISH methodologies.

IMAGING AND REPORTABLE RANGE




Figure 5: Representative images taken for validation samples. Image on right is from NCI-H827 cell line positive for 8x amplification of PDGFR α amplification and the image on the left is negative for both genes of interest.

Reportable ranges: A cohort of 23 known negative samples was used to calculate assay cut- off for PDGFR α copy number at 3.65 copies. Samples that have >6 PDGFR α average copy number or a test over reference ratio of ≥ 2.0 are classified positive. Samples with a PDGFR α copy number between 4 and 6 are classified as equivocal /low positive if the ratio is <2.0. RET gene rearrangements usually occur in high frequencies and the calculated cut off for classical 1R1G1F pattern was 4.0%. All other rearrangement patterns will be evaluated at greater than 5.0% based on significance.

REPEATABILITY AND REPRODUCIBILITY

Sample ID	RET Break-apart %/ Avg PDGFR α copy No. Repeat 1/ Day 1	RET Break-apart %/ Avg PDGFR α copy No. Repeat 2/ Day 1	RET Break-apart %/ Avg PDGFR α copy No. Repeat 3/ Day 1	RET Break-apart %/ Avg PDGFR α copy No. Repeat 4/ Day 2	RET Break-apart %/ Avg PDGFR α copy number Repeat 5/ Day 5	Standard deviation RET Break-apart/ Avg PDGFR α copy No.	% Confidence Interval
GM19204	2% / 2.26 copies	2% / 2.30 copies	2% / 2.35 copies	0% / 2.20 copies	2% / 2.32 copies	0.014/ 0.06	0% / 2.54%
LC2Ad	68% / 3.16 copies	60% / 3.20 copies	62% / 3.24 copies	54% / 3.1 copies	60% / 3.5 copies	0.05/ 0.15	8.25%/ 4.76%
Sample I	2% / 2.64 copies	2%/ 2.68 copies	2% / 2.81 copies	2% / 2.70 copies	2% / 2.75 copies	0.07/ 0	0% / 2.42%

Figure 4: Above are the results of the assay repeatability and reproducibility criteria. A set of three samples were selected from the validation dataset and repeated three times on the same day and twice more on different days. The results were combined and standard deviation calculated to derive the percent co- efficient of variation. The target criteria for assay acceptance is >95% concordance with reference results and percent co-efficient of variation of <25%. Both criteria have been met above successfully. The entire validation was performed on at least 20 normal/ negative samples from Non- Small Cell Lung Carcinoma cell lines or Human Bio specimens in order to derive cut- offs for PDGFR α amplification and RET gene rearrangement.

REFERENCES

- Caria P, Frau DV, Dettori T, et al. Optimizing detection of RET and PPAR γ rearrangements in thyroid neoplastic cells using a home-brew tetracolor probe. *Cancer Cytopathol.* 2014;122(5):377-85.
- Caria P, Dettori T, Frau DV, et al. Assessing RET/PTC in thyroid nodule fine-needle aspirates: the FISH point of view. *Endocr Relat Cancer.* 2013;20(4):527-36.
- Fu S, Liang Y, Lin YB, et al. The Frequency and Clinical Implication of ROS1 and RET Rearrangements in Resected Stage IIIA-N2 Non-Small Cell Lung Cancer Patients. *PLoS ONE.* 2015;10(4):e0124354.
- RET-rearranged non-small-cell lung carcinoma: a clinicopathological and molecular analysis. *British Journal of Cancer.* 2014;110(6):1571.
- Wang R, Hu H, Pan Y, et al. RET fusions define a unique molecular and clinicopathologic subtype of non-small-cell lung cancer. *J Clin Oncol.* 2012;30(35):4352-9.
- Shaw AT, Ou SH, Bang YJ, et al. Crizotinib in ROS1-rearranged non-small-cell lung cancer. *N Engl J Med.* 2014;371(21):1963-71.
- Ramos, A. H., Dutt, A., Mermel, C., Perner, S., Cho, J., Lafargue, C. J., . . . Meyerson, M. (2009, November 01). Amplification of chromosomal segment 4q12 in non-small cell lung cancer. *Cancer Biology & Therapy Cbt, 8(21), 2042-2050.* doi:10.4161/cbt.8.21.976